

XDAQ AIO

Updated Gen2 Design

Optimized Hardware, Firmware, API Fast, Low Latency Interface

Better Signal, All Electrically Isolated Ports Powerful, Yet Portable

Neuropixels Out-of-box



State-of-the Art Neuropixels Support - Seamlessly integrated support for both Neuropixels 1.0 and 2.0, enabling high-density recordings of up to 3,072 channels per unit.

High-Channel Passive Probe Compatability - Capture EEG, ECoG, LFP, and action potentials through a wide range of silicon probes or microwires - supporting up to 256 channels per cable and 1,024 channels per box.

Built-in Electrical Stimulation - Delivers up to 64 channels per cable and 128 channels per unit - fully integrated for streamlined, high-fidelity neural stimulation

Integrated Optogenetics - Built-in 4-channel high-power driver for precise, reliable optical stimulation - no external hardware required.

Robust I/O, HARP-Compatible, and Scalable - Each XDAQ features 32 digital I/O pairs and 8 ADC/DAC pairs. HARP-compatible as a master, it integrates with existing HARP systems and links effortlessly with other XDAQ.

Guaranteed Synchronized Behavior Video Capture - AlO's onboard video engine precisely aligns neural and video data by embedding neural timestamps into video frame metadata. Supports multiple high-speed cameras (>120 FPS) for multi-angle monitoring or standard cameras for multi-animal monitoring.



EIB & Adapter

Off-the-shelf adapters for seamless integration with popular opensource platforms and commercial probes.



Cable

Optimized, cost-effective cabling solutions designed to reduce the long-term cost of in vivo experiments across Neuropixels and conventional probe technologies.



Validation Tool

A suite of tools including brainwave-compliant signal generators, stimulation pulse validators, and impedance test loads to verify system integrity and ensure experimental reproducibility.



IO and Headstage Expansion

Modular expansion options to increase digital I/O channels, support additional independent headstages (for large scale electrophysiology monitoring), or enable multiple XDAQ systems to operate in synchronized tandem.

New: HARP Adapter now available



Commutator

Low-noise, cost-effective, and stress-tested motorized commutators compatible with both Neuropixels and Intan ecosystems.


A non-motorized version is also available for the Intan platform.



MEA Interface

Cost-effective MEA accessory suite featuring a moisture-resistant design and microscopy access window for seamless XDAQ integration. The customizable headstage supports both stimulation and recording, and is compatible with a variety of MEA configurations.

For more details, you can scan this code to view the extended flyer:





Fullspec Passive Probe Recording (XR Mode)

Capture EEG, ECoG, LFP, and action potentials using a wide range of silicon probes or microwires.

Recording capability (when paired with XR(RHD) Headstages:

16-bit, 30kS/s/ch

±5 mV input range

2.4 μVrms input-referred noise
 15 μV min detectable signal

On-chip DSP high pass filter 100-20 kHz (upper cut off filter) 0.1-500 Hz (lower cut off filter

in-situ impedance measurement

Fullspec Passive Probe Recording (XR Mode)

Deliver constant-current stimulation and record from the same amplifiers. XDAQ ONE and AIO systems enable switchable VStim compliance levels, providing up to 40% additional peak current.

Stimulation Capability (when paired with XSR(RHS) Headstages:

Switchable rec/stim amplifier

Constant current driver

10nA-2.55 mA 14V compliance

Flexible stimuli, 33 µs step

Biphasic Triphasic

Burst

Specifications	XR Mode	XSR Mode	NP M	ode
Neural Acquisition			NP 1.0	NP 2.0
Channel Count	1024ch	128ch	1536ch	3072ch
Sampling Rate	1-30kHz		AP: 30 kHz LFP: 2.5kHz	30kHz
A/D Resolution	16-	-bit	10-bit	12-bit
Input Range	±5	mV	±5 mV	10 mV(ppk)
Input-referred Noise (RMS)	2.4 μV	typical	AP: 5.9 μV typical LFP: 9.2 μV typical	' ''
Hardware Filter		ff: 100 - 20 kHz ff: 0.1 - 500 Hz	AP: 0.3-10 kHz LFP: 0.5-500 Hz	0.5-10 kHz

Electrical Stimulation	XSR Mode	NP Mod	de		
Channel Count	128 ch				
Stim Magnitude, Type	10n to 2.55mA, constant current				
Stim Pulse	Biphasic, Triphasic, Burst	Biphasic, Triphasic			
Stim time resolution	33 µs				
Compliance	±7V, -4 to +10V, -10 to +4V				
Synchonized GPIO					
Digital In (2.2 - 5.5V Logic High	8ch onboard	32ch expandable			
Digital Out (BNC 3.3, 5V configural (DB25: 3.3V)	le) 8ch onboard 32ch expandable				
ADC (16bit, ± 10 V)	2ch onboard 8ch expandable				
DAC (16bit, ± 10 V)	2ch onboard 8ch expandable				
Connector and Port Headstage port (electrically iso Peripheral & System GPIO Video System Optogenetics LED Driver Audio Monitoring IO Expansion Communication Power Ground Terminals	8x E	IDMI Type A BNC Connectors, 1X D 4X USB Type A, U 1X D-Sub9 Co 1X 3.5mm 2X HDMI Ty 2X Thunderbolt 3 (USI 1X DC Jack(2.5mm, 5) 1X Chasis GND, 1X	VC camera nnector Jack pe D B-C Connector) 5525) 20V 10A		
Dimension (mm)	227x190x94				
Weight (g)		3045			
Software	XR Mode Radiens, OpenEphys XDAQ™ plugin, XDAQ-RHX	XSR Mode Radiens, XDAQ-RHX	NP Mode OpenEphys XDAQ™ Neuropixels plugin		